Device & consumable description

Juvapen (Toxin)

The Juvapen (Toxin) is a cordless motorized injection system powered by a lithium battery. The Juvapen (Toxin) is a non-sterile device intended to assist certified practitioners in injecting Botulinum toxins. After attaching the syringe holder (Juvalink (Toxin)) onto the device, the practitioner simply inserts the dedicated filled 1mL syringe onto the Juvalink (Toxin). Once the Juvapen (Toxin) is fully assembled the practitioner presses on the top part of lever or the syringe holder for an accurate dosage delivery in drops or continuous flow mode. The Juvapen (Toxin) is only intended for Botulinum toxin injections and should be used only on adult patients.

Syringe holder / Juvalink (Toxin)

The Juvalink (Toxin) has been designed to improve syringe handling ergonomics, allowing the user to have greater freedom to operate during injections. The Juvalink (Toxin) provides low force accurate activation with the lever or strong grip activation directly pressing on the top of the syringe holder. The Juvalink (Toxin) is only compatible with the 1 mL syringe that is provided. The Juvalink (Toxin) is compatible exclusively with Juvapen (Toxin).

Advantages of the Juvapen (Toxin)

For the patient:

- Significantly less pain.
- Highly precise injections.
- Reduced side effects.

For the practitioner:

- Easy to use & intuitive system.
- Optimized control of the procedure: depth of the injection, homogeneous delivery and precise dosage.
- Accuracy in the dosages.

Precautions and recommendations

Precautions for use: please read the following instructions carefully and follow the proper safety measures in order to perform a successful injection.

- Only a qualified and experienced medical practitioner should use the Juvapen (Toxin)
- The use of the Juvapen (Toxin) and the results thereof are the exclusive responsibility of the practitioner.
- The practitioner is the only one who can evaluate the indications and contraindications relative to the technique to employ, the product to be administered and the appropriate needle to be used.
- In all cases, the medical practitioner shall proceed only after having fully familiarized himself / herself with both the product and its educational material, including their instructions for use.
- Only 1ml syringes that are provided with are recommended by the manufacturer.
- The Juvapen (Toxin) syringe and its holder are for a single patient and single session use. Discard immediately after use into a sharps waste container. Do not use if the package is opened or damaged.

Adverse reactions and contraindications

It is the practitioner's full responsibility to read and follow the Instructions for use in order to avoid any adverse reactions or complaints. Possible side effects and contraindications can result from the use of products in conjunction with the Juvapen (Toxin) (e.g. Botulinum toxin, needle, cannula, etc.). It is the practitioner's full responsibility to check and follow the instructions for use of any products used with the Juvapen (Toxin).

General recommendations

- Keep the Juvapen (Toxin) out of reach of children.
- The equipment cannot be used where there is an inflammable mixture of air, oxygen or nitrogen.
- Store the Juvapen (Toxin) away from heat and direct exposure to sun.
- Keep the Juvapen (Toxin) away from water and chemical products.
- Always check whether the device is functioning properly before use.
- If the Juvapen (Toxin) is damaged, it must be sent back to Juvaplus SA (please refer to the Warranty).
- If the On/Off light turns orange, change the battery.
- The Juvapen (Toxin) turns off automatically, when not in use, after 2 minutes.
- The Juvalinks (Toxin) are designed for use with the Juvapen (Toxin) exclusively and therefore cannot be used with any other third party device.
- Do not use other batteries other than the ones supplied by Juvaplus SA. Damage of the Juvapen (Toxin) caused by third party batteries voids the warranty.
- Do not sterilize any of the components: Juvapen (Toxin), syringe holder and battery.
- Wearing single-use gloves is highly recommended for the use of the equipment.
- No alteration of this equipment is allowed.
- Dispose of the battery according to your local environmental regulations.
- Do not dispose of the Juvapen (Toxin), but rather, return it to Juvaplus SA.

- Do not use the Juvapen (Toxin) near electronic equipment which emits electromagnetic fields.
- Only use syringes provided by Juvaplus SA.

Accessories

Battery cap Lithium 6V /170 mAh

The Juvapen (Toxin) must be used with a certified battery IEC60086-4 only Juvalink (Toxin)

Juvalink (Toxin): Plastic PC Makrolon 2458

Cleaning

After each use, the Juvapen (Toxin) must be cleaned with an adequate cleaning agent before re-using it. Do not immerse the device or its components in water. We recommend using an appropriate cleaning product in accordance with your local regulation.

Handling conditions

It is recommended to store and transport the Juvapen (Toxin) in its original packaging.


	OPERATING 0 - 2000 m	TRANSPORT & STORAGE
Temperature	25°C 77°F 	40°C 104°F 
Humidity	60% RH at 86°F / 30°C non-condensing	

Assembly and use

Please refer to the assembly instructions drawings attached.

1. The battery is inserted by screwing it into the Juvapen (Toxin) (do not touch the keyboard while screwing).
2. Unpack the syringe and fill it with toxin in the desirable volume.
3. Break the plunger rod at the chosen volume of toxin.
4. Remove all air bubbles that might have accumulated in the syringe.
5. Connect the Juvalink (Toxin) with the Juvapen (Toxin) until you hear the Click that indicates the proper attachment with the device.
6. To attach the syringe to the Juvalink (Toxin), rotate the syringe onto the bayonet mount. It is properly attached when it is aligned with the device and you feel a firm click.
7. Turn on the device by pressing the button (bottom green light indicates the device is switched on)
8. Attach the needle to the syringe and remove its protection.
9. Select the continuous flow mode, indicated by an orange triangle and orange light.
10. Remove the air from the syringe and prime the needle by pressing on the activation lever, keep on pressing until it removes all the air.
11. Select the desirable dosage by choosing the correct volume (blue led indicates the selected dosage in milliliter).
12. The injection can be performed by pressing the Juvalink (Toxin) lever and the Juvapen (Toxin) will execute the dosage option.
13. Replace the needle protection.
14. Rotate the syringe from the syringe holder in order to remove it and then discard the syringe and the needle.
15. Select the rewind mode (indicated by an upside down orange arrow and orange light) – If you have completely injected the contents of the syringe, Juvapen will automatically go into Rewind mode.
16. Press Juvalink (Toxin) lever once to rewind the lead screw to its original position.

The Device Control Buttons and Control Lights

Warranty

The warranty period for the Juvapen (Toxin) is 24 months from the date of purchase.

Status light color	Continuous or Flashing	Indication	Recommended action
Green	Continuous	Ready to use	Normal use of the device
Orange	Continuous	Warning: End of battery life	Consider changing battery
Orange	Flashing	Piston to its end-point <i>Occlusion of the needle</i>	Juvapen will automatically go into Rewind mode. Press Juvalink lever once to return the lead screw to its original position. <i>Change the needle</i>
Red	Flashing	Defective sensor	Contact Juvaplus
Red	Continuous	End of battery life	Change the battery

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
Juvaplus Europe
15 rue des Halles,
75001 Paris
FRANCE
info@juvaplus.com

EMC information

JUVAPEN complies with the standard IEC 60601-1-2: 2014, Class B/Group 1 in emission.
Medical electrical equipment should be used with precautions according to EMC notices disclosed in this manual.
JUVAPEEN is intended for use in the electromagnetic environment specified below. The customer or user of JUVAPEN should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
Conducted emissions CISPR 11: 2009 / AMD1: 2010	Not applicable	JUVAPEEN uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. JUVAPEEN is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Radiated emissions CISPR 11: 2009 / AMD1: 2010	Class B / Group 1	
Harmonic emissions IEC 61000-3-2:2018	Not applicable	
Voltage fluctuations flicker emissions IEC 61000-3-3: 2013/AMD1:2017	Not applicable	
Devices with motors or switching circuits CISPR 14-1	Not applicable	

Immunity test	IEC60601 test level	Verdict	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2: 2008	+/- 8 kV contact discharge +/- 15 kV air discharge	Pass	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4: 2012	+/- 2 kV AC Mains +/- 1 kV I/O Ports 100 kHz PRR	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. JUVAPEEN is a self-powered (non-rechargeable battery) handy device that is not connected to the public low-voltage power supply network.
Surge – AC Mains IEC 61000-4-5: 2014 AMD1: 2017	+/- 2 kV (max) 0,90, 180 & 270 ° phase angles	Not applicable	
Surges – 12 VDC Power ISO 7637-2	600 V	Not applicable	
Voltage dips & short interruptions immunity IEC 61000-4-11: 2004 AMD1:2017	UT=0%, 0,5 cycle (0, 45, 90, 135, 180, 225, 270 & 315°) UT=0%; 1 cycle – UT=70% 25/30 cycles (0°) UT=0%; 250/300 cycles	Not applicable	
Conducted Immunity IEC 61000-4-6:2013	3V (0.15-80MHz) 6V (ISM bands)	Not Applicable	
Magnetic Immunity IEC 61000-4-8: 2009	30 A/m 50 and 60 Hz	Pass	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC60601 test level	Verdict
Radiated RF IEC IEC 61000-4-3: 2006 AMD1:2007/AMD2:2010	10V/m 80 MHz to 2,7GHz	Pass
Proximity Field from Wireless Transmitters IEC 61000-4-3: 2006 AMD1:2007/AMD2:2010	9V/m to 28 V/m 15 specific frequencies	Pass
<p>Electromagnetic environment – guidance Portable and mobile RF communications equipment should be used no closer to any part of the JUVAPEN, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance $d = 1.17 \sqrt{P} / 0.15 \text{ MHz to } 80 \text{ MHz}$ - $d = 1.17 \sqrt{P} / 80 \text{ MHz to } 800 \text{ MHz}$ - $d = 2.33 \sqrt{P} / 800 \text{ MHz to } 2,7 \text{ GHz}$, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>		

Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz à 80 MHz $d = 1.17 \sqrt{P}$	80 MHz à 800 MHz $d = 1.17 \sqrt{P}$	800 MHz à 2,7GHz $d = 2.33 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.33 m
10	3.7 m	3.7 m	7.37 m
100	11.70 m	11.70 m	23.30 m