

# DENTAPEN<sup>®</sup>




**DENTAL ANESTHESIA**  
Painless Effortless Cordless

# FOR YOUR PATIENTS

“The number 1 criteria in choosing a dentist is:  
A dentist who can provide painless anesthesia.”

*How dentists are judged by patients, St Georges J, Dent Today, 2004.*

Pressure  
(N/mm<sup>2</sup>)



Controlled  
Injection

=

Increased  
Success Rate

+

Reduced  
side-effects

+

Minimized  
numbing

# DENTAPEN®

ROTATING BARREL  
for easy bevel placement



AUTOCLAVABLE  
low operating costs



BATTERY  
up to 100 patients

2-YEARS WARRANTY


Complete Dentapen® Kit:

Metric hub ..... Ref. DK01.2/M

Imperial hub ..... Ref. DK01.2/I

Optional highly recommended:

Protective sleeve x 250 ..... PS01



- 1 Dentapen® x1
- 2 Finger grip x5
- 3 Finger grip for nerve block x2
- 4 Cartridge holders x7  
1.7ml or 2.2ml  
metric or imperial
- 5 Battery x3

Includes:  
protective sleeves

# Designed for EFFICIENCY

- ✓ Autoclavable accessories
- ✓ Universal: all techniques
- ✓ Two-way activation and grip: syringe or pen like
- ✓ Battery powered
- ✓ Compact and light: only 50 grams
- ✓ Active aspiration



## Easy Set-up



- ☀ Slow 90 s/mL
- ☀ Medium 60 s/mL
- ☀ Fast 30 s/mL



- ☀ Intraligamentary mode
- ☀ Ramp-up mode = anesthesia flow increases gradually for totally painless procedure



Rewind back to initial position



# FOR YOU

## Cutting-edge C-CLAD\* technology for advanced and painless dental anesthesia

- ✓ Controlled initial force and constant flow for smoother injections:
  - Increases anesthesia success rate and saves time
- ✓ Computerized pressure feedback sensor:
  - Adjusts the injection speed depending on pressure of the tissue
- ✓ Other benefits:
  - Allows challenging injections: intraligamentary
  - Performs aspirations for a safe procedure
  - Decreases risk of work related injuries: tendinitis, arthritis
  - Manages patient anxiety



## LESS



PAIN



ANXIETY



NUMBING

## MORE



CHILDREN  
FRIENDLY



SATISFIED  
PATIENTS



REFERRAL

\*computer-controlled local anesthesia device



[www.juvaplus.com](http://www.juvaplus.com)

 @dentapen\_injector

 Juvaplus

 Dentapen

R O B O T I C   S Y R I N G E S

J U V A +

Designed and manufactured in Switzerland  
Chemin Louis-Hubert 3, CH-1213 Petit-Lancy  
[info@juvaplus.com](mailto:info@juvaplus.com) - Tel. +41 76 286 34 37